

LIMNIC ERUPTIONS STATION #3A

KWLA Question: Why did the water sample from the bottom of the lake bubble when brought to the surface?

Investigation: How does the environment at the bottom of the lake vary from the top in terms of pressure?

Materials:

- Vernier LabQuest⁹
- Vernier Gas Pressure sensor⁹
- Fluorescent bulb cover (1.2 m) glued to a PVC cap with silicone caulk
- Flexible tubing (1.5 m) containing a weight (fishing sinker) attached to one end with a paper clip
- Bin to catch water overflow

Procedure:

1. Place tube in plastic bin to catch any spilled water.
2. Fill tube with water to the 0 cm mark.
3. Attach pressure sensor to flexible tubing.
4. Using LabQuest to collect gas pressure data every 10 cm.
 - a. Tap Mode.
 - b. Tap the arrow beside the top bar that reads time graph.
 - c. Select Events with Entry.
 - d. In name box, type “depth” and in units box type “cm”.
 - e. Tap OK.
 - f. Press the START button.
5. Place the opening of the flexible tubing at the 0 cm mark directly above the water and allow the reading to stabilize. **ALWAYS KEEP THE PRESSURE SENSOR ABOVE THE TUBE SO THAT WATER DOESN'T GET INTO THE SENSOR.**
6. Collect data.
 - a. Tap Keep.
 - b. Type depth of “0” cm.
 - c. Tap OK.
7. Drop the end of the tubing to the 10 cm mark and allow the reading to stabilize. Tap Keep, type in a depth of “10” cm, and tap OK.
8. Repeat step #7 for every 10 cm.
9. Sketch graph.

Outcome: (Analyze in your Limnologist's Journal.)

*See Sources page for footnote references.