# HYDROFLUORIC ACID

#### **Hazard Description**

Hydrogen fluoride (HF) is used at the University of Waterloo in various forms. Hydrofluoric acid comes in liquid form in various concentrations. It is clear, colourless, and resembles water. Anhydrous hydrogen fluoride is a gas at room pressure (1 atm) and temperature (20°C). Regardless of its physical state or concentration, it is a highly toxic and corrosive substance with an irritating, pungent odour.

### Prior to Working with Hydrogen Fluoride (HF)

- Complete a research-specific laboratory risk assessment.
- Ensure that a standard operating procedure (SOP) is created, posted, trained and approved by the supervisor on the process.
- Only use in a fume hood.
- Do not use HF when working alone or outside regular business hours (8 am 4:30 pm).
- Do not work with HF until you have received practical training on all the of the specific work procedures of your laboratory (or workspace) that use HF. This should be performed by someone knowledgeable in the use and hazards of HF. This includes your supervisor, a lab technician, or senior graduate student.
- Determine if any suitable substitute can be found. If no substitute is available and HF is required for the work being undertaken, the Safety Office requires the lab purchase 2.5% calcium gluconate gel (contact ChemStores). Note this has an expiry date, and all users need to have this ready for immediate application before beginning work and while storing HF.
- Do not purchase HF in larger than 2L. Discard any HF that has not been used for the last 6 months.
- Post a sign where HF will be used in the work area and say "Hydrogen Fluoride" not "HF."
- Ensure that all required materials are available in case of spill, emergency, and for waste collection.
- Ensure an HF-specific spill kit is available which includes everything in a basic spill kit plus:
  - HF-specific neutralizer
    - Not all HF neutralizers are made the same. If you would like to purchase a
      different neutralizer, please have it reviewed by <a href="mailto:chemsafety@uwaterloo.ca">chemsafety@uwaterloo.ca</a> in
      advance.
  - Absorbent, appropriate for strong acids

#### Spill response

- 1. Clear the area and call the Spills Team.
- 2. Lay absorbents on the area to prevent vapours.
- 3. Using gloves, pick up the saturated absorbent, and continue until the liquid is absorbed.
- 4. Using tweezers, pick up any remaining plastic or equipment pieces involved in the spill.
- 5. Spray the area with the HF neutralizers, and briefly let it sit. Check the area with pH paper.
- 6. Place absorbent on the area to pick up the neutralizer.
- 7. Retest with pH paper to ensure the spill is no longer acidic.
- 8. Dispose of as Hydrofluoric Acid contaminated waste.
- 9. If it caused damage to the floor, contact custodial to have the issues addressed.
- 10. You can damp paper towels to remove any sticky residue or neutralizer and the water can be disposed of in the drain.

#### Handling

- Before using any HF, inspect the container to ensure it is not defective or damaged.
- All work materials, including syringes, should be made of compatible material for HF. HF dissolves most metals, natural rubber, concrete, glass, fiberglass, ceramics and glazes. HF does not attack metallic lead and platinum, polyethylene, polypropylene, Teflon, Plexiglas (i.e., acrylic), and wax.
- Additional PPE is required when working with HF:
  - Gloves: either double gloved Nitrile gloves, PVC, or neoprene gloves, or a combination of the two (one layer of nitrile gloves is not sufficient for working with HF)
  - Long pants, long sleeves, lab coat and closed-toe shoes
  - Goggles
  - Face shield
  - Acid -resistant apron

## **Emergency Procedures**

Always review SDS of purchased product for manufacture specific recommendations. Look at SDS for other modes of exposure.

Depending on the concentration of exposure, the effect of HF can be delayed for up to 24 hours so medical attention should always be received after HF exposure.

Concentration	Time to Onset of Symptoms	Description of Surface Reaction	System Toxicity
> 50 %	Immediate	Immediate burns, rapid destruction of tissue, severe pain.	Free fluoride ions aggressively bind to cations such as calcium, magnesium, and potassium causing:
20 % - 50 %	1 to 8 hour delay	Burns and destruction of tissue with pain	<ul> <li>Interference with nerve impulses, muscle contraction and relaxation</li> <li>Electrolyte imbalance leading to an</li> </ul>
< 20 %	Up to 24 hour delay	Redness, mild to severe pain, or potentially no surface reaction.	<ul><li>irregular heartbeat or heart attack</li><li>Lowering of blood pH (metabolic acidosis)</li></ul>

When exposed to hydrogen fluoride gas, immediate effects include:

- Upper respiratory tract irritation (coughing, choking)
- Difficulty breathing, chest tightness

Contact by acid solution potential immediate effects:

- Burns to the eyes, opacity of cornea
- Burns to the skin

In many cases, exposure to HF solutions with low concentrations can lead to delayed reactions.

## No concentration is safe to handle without appropriate precautions.

	Contacts		
	Emergency: 911		
	UW Special Constables: 519-888-4911 or ext. 22222		
	Poison Control: 1-800-268-9017		
Whenever 911	is called, if possible, UW Special Constables should also be informed to make them aware of the		
emergency on o	campus and allow them to support as needed. Ask them to meet the paramedics and direct them to		
the incident loca	ation.		
Inhalation	Remove individual from contaminated area		
	Call 911 for transport to hospital		
	Corrosive substances may cause severe lung damage if inhaled		
	Perform CPR and artificial respiration if necessary		
	Provide responder with SDS		
Skin Contact	Call 911 for transport to hospital and inform them of HF exposure		
	Remove contaminated clothing and quickly but gently wipe material off skin		
	<ul> <li>Flush with water and apply calcium gluconate (can be applied during shower)</li> </ul>		
	Apply calcium gluconate every 15 minutes and can continue while people transported to the		
	hospital		
	Provide responder with SDS		
Eye Contact	Call 911 for transport to hospital		
	Flush eyes using eyewash station for a minimum of 15 minutes		
	Do not apply calcium gluconate to the eyes		
	Provide responder with SDS		

### Storage

Storage Group F – Inorganic Acids (with further segregation from other inorganic acids)

#### **Hazardous Waste**

- Must be disposed of as per the <u>University's Hazardous Waste Standard</u>.
- Inorganic acids should not be mixed prior to being sent to the waste facility. Individually segregate.
- Label and make a waste Hydrofluoric Acid container that will only be used for hydrofluoric acid. The waste hydrofluoric acid container should only be made of polyethylene or Teflon.
- If bottles are degraded, contact <u>esf@uwaterloo.ca</u> for disposal instructions.